
















































































































































Straumann® Dental Implant System System Overview TISSUE LEVEL PROSTHETICS 1.0 גרסה 2019	NNC Ø 3.5 mm  CrossFit®	RN Ø 4.8 mm  synOcta®	WN Ø 6.5 mm  synOcta®	NNC Ø 3.5 mm 	RN Ø 4.8 mm 	WN Ø 6.5 mm 			
Impression Posts Impression Caps Implant Analogs	Open-tray Impression  048.128 Closed-tray Impression  048.134 Analog  048.127	Implant Level Open-tray Impression  048.090/048.010 Closed-tray Impression  048.862 Analog  048.124 Abutment Level Open-tray Impression  048.866 Closed-tray Impression  048.864 Analog  048.108	Implant Level Open-tray Impression  048.091 Closed-tray Impression  048.870 Analog  048.171 Abutment Level Open-tray Impression  048.869 Closed-tray Impression  048.871 Analog  048.872						
Digital Impression Posts Scanbodies and Repositionable Analogs	 048.173	 048.129  048.168	 048.172  048.169						
Temporary Restorations Protective Caps	 •048.501 Crown incl. screw 048.313  ••048.502 Bridge incl. screw 048.313  048.699 Crown incl. screw 048.314	 •048.715 Crown incl. screw 049.181  ••048.716 Bridge incl. screw 049.181  •048.664 incl. screw 048.356	 •048.717 Crown incl. screw 049.181  ••048.718 Bridge incl. screw 049.181  •048.242 incl. screw 048.356						
Prosthetic Restoration		Screw-retained	Screw-retained	Cement-retained	Screw-retained	Cement-retained	Locator®	Locator®	Locator®
Titanium Abutments	Cement-retained  048.598 incl. screw 048.314  048.621 incl. screw 048.314 Angled  048.622 incl. screw 048.314 Screw-retained Variobase® for Crown  AH 3.5 mm incl. screw 048.313  AH 5.5 mm incl. screw 048.313 Variobase® for Crown AS  AH 3.8 mm incl. screw 048.893  AH 5.8 mm incl. screw 048.893 Variobase® for bridge/bar  AH 3.5 mm incl. screw 025.2926 Variobase® C  incl. screw 048.313	Screw-retained  synOcta® 1.5 048.601 Variobase® for Crown  AH 3.5 mm incl. screw 048.356  AH 5.5 mm incl. screw 048.356 Variobase® for Crown AS  AH 4.0 mm incl. screw 048.894  AH 6.0 mm incl. screw 048.894 Variobase® for bridge/bar  AH 3.5 mm incl. screw 048.356 Variobase® C  incl. screw 022.0045	Cement-retained  022.0034 TAN CAD/CAM ¹⁾ 2)	Screw-retained  synOcta® 1.5 048.603	Screw-retained  synOcta® Angled 048.612/613/617/618 H 5.7 mm 048.610/611/615/616 H 6.7 mm  synOcta® Cementable incl. screw Screw-retained Variobase® for Crown  AH 3.5 mm incl. screw 048.356  AH 5.5 mm incl. screw 048.356 Variobase® for Crown AS  AH 4.5 mm incl. screw 048.894  AH 6.5 mm incl. screw 048.894 Variobase® for bridge/bar  AH 3.5 mm incl. screw 048.356 Variobase® C  incl. screw 022.0045	Cement-retained  022.0036 TAN CAD/CAM ¹⁾ 2)	Locator®  H1 mm 048.581  H2 mm 048.582  H3 mm 048.583  H4 mm 048.584  H5 mm 048.585  H6 mm 048.586 Up to 20° Divergence (Pack of 10)  8519-10	Locator®  H1 mm 048.175  H2 mm 048.176  H3 mm 048.177  H4 mm 048.178  H5 mm 048.179  H6 mm 048.180 Up to 20° Divergence (Pack of 10)  8519-10	Locator®  H1 mm 048.183  H2 mm 048.184  H3 mm 048.185  H4 mm 048.186 H5 mm 048.187 Up to 20° Divergence (Pack of 10)  8519-10
Plastic Copings	 •048.256 ••048.257  048.258 Burn-out Coping for Variobase® for Crown  AH 3.5 mm 048.267*  AH 5.5 mm 023.0014*  048.896 Burn-out Coping for Variobase® for Crown AS  AH 3.5 mm 048.380	 •048.229 Crown  ••048.227 Bridge excl. screw 048.350 כניע ללא בורג Burn-out Coping for Variobase® for Crown  AH 3.5 mm 048.268*  AH 5.5 mm 023.0015*  048.897 Burn-out Coping for Variobase® for Crown AS  AH 3.5 mm 048.268*	 048.670  048.676/4  ••048.662 Bridge  •048.663 Crown	 •048.667 excl. screw 048.350 כניע ללא בורג Burn-out Coping for Variobase® for Crown  AH 3.5 mm 048.269*  AH 5.5 mm 023.0016*  048.898 Burn-out coping for Variobase® for Crown AS  AH 3.5 mm 048.382	 048.678  •048.243  ••048.244	 8540-10 Up to 40° Divergence (Pack of 10)	 8540-10 Up to 40° Divergence (Pack of 10)	 8540-10 Up to 40° Divergence (Pack of 10)	
Gold Copings Gold Abutments	 •048.592  ••048.593 incl. screw 048.313 Gold Abutment	 •048.633 Crown  ••048.632 Bridge excl. screw 048.350 כניע ללא בורג synOcta® Gold Copings	 •048.642  ••048.802 incl. screw 048.356 synOcta® Gold Abutment	 •048.639  ••048.638 excl. screw 048.350 synOcta® gold coping	 •048.644  ••048.803 incl. screw 048.356 synOcta® gold coping	 048.356  022.045  Basal Screw AS 048.894  048.356			
Auxiliary Parts Screws	 048.314 Cementable  048.313 VB + Gold  025.2926 VB Bridge  048.893 AS Basal Screw AS	 048.350  048.356 VB  022.0045  048.894 Basal Screw AS	 048.350  048.356	 048.350  048.356	 048.356 VB Crown  022.045 VB C  048.894 Basal Screw AS	 048.356		נפתח פרוטסי  046.400 15 mm  046.401 21 mm  046.402 27 mm  046.119  046.049	