

### RB/WB prosthetics – fit on RB and WB implants (∅ 3.5 mm - ∅ 6.5 mm)

• crowns with slim emergence profiles • bridges

### Straumann® BLX Implant TorcFit™ Connection

#### RB / WB ∅ 3.5 - ∅ 6.5

**Implant abutment compatibility**

RB/WB and WB  
RB/WB Variobase® for Crown

compatible with

RB implants ∅ 3.5 - 4.5 mm  
WB implants ∅ 5.0 - 6.5 mm

**Warning:** No WB abutments on RB implants!

BLX RB (Regular Base) ∅ 3.5 mm

BLX RB (Regular Base) ∅ 3.75 mm

BLX RB (Regular Base) ∅ 4.0 mm

BLX RB (Regular Base) ∅ 4.5 mm

BLX WB (Wide Base) ∅ 5.0 mm

BLX WB (Wide Base) ∅ 5.5 mm

BLX WB (Wide Base) ∅ 6.5 mm

#### RB / WB ∅ 3.5 - ∅ 6.5

**Implant abutment compatibility**

RB/WB and WB  
RB/WB Variobase® for Crown

compatible with

RB implants ∅ 3.5 - 4.5 mm  
WB implants ∅ 5.0 - 6.5 mm

**Warning:** No WB abutments on RB implants!

#### WB ∅ 5.5 - ∅ 6.5

+WB+ WB Variobase® for Crown

compatible with

#### Temporary Abutments for Crown, AH 11 mm

		∅ 3.8 mm (RB/WB)	∅ 4.5 mm (RB/WB)	∅ 6.0 mm (RB/WB)
GH	1.5 mm	064.4362	064.4372	—
	2.5 mm	064.4363	064.4373	064.4382
	3.5 mm	064.4364	064.4374	064.4383
<b>Incl. Screw</b>		065.0036		

#### Temporary Abutments for Bridge, AH 5.0 mm

		∅ 4.5 mm (RB/WB)
GH	1.5 mm	064.4352
<b>Incl. Screw</b>		065.0036

#### Open-tray Transfer Closed-tray Transfer Implant Analog Metal Scanbody Scanpost for Sirona® Cerec Digital Analog

16mm 065.0031

24mm 065.0033

13mm 065.4310

065.0021

065.0103

065.0038

065.0023

#### SRA Closed-tray Impression for Screw-retained Abutments

		∅ 4.6 mm
Crown		025.2246
Bridge/Bar		025.0014

#### SRA Open-tray Impression for Screw-retained Abutments

		∅ 4.6 mm
Crown		025.2244
Bridge/Bar		025.0012

**SRA Analog**

Abutment Analog		023.4756
-----------------	--	----------

**Straumann® Metal Scanbody for Screw-retained Abutment**

	∅ 4.6 mm	025.0081
--	----------	----------

**SRA Analog - Digital**


Repositionable Analog		025.0008
-----------------------	--	----------

#### Temporary Abutments for Crown, AH 11 mm

		∅ 5.5 mm (WB)
GH	0.75 mm	064.4391
	1.5 mm	064.4392
<b>Incl. Screw</b>		065.0036

#### Temporary Abutments for Bridge, AH 5.0 mm

		∅ 4.5 mm (RB/WB)
GH	1.5 mm	064.4352
<b>Incl. Screw</b>		065.0036



**Variobase® for Crown, AH 5.5 mm**

		∅ 3.8 mm (RB/WB)	∅ 4.5 mm (RB/WB)
GH	1.5 mm	062.4934	062.4944
	2.5 mm	062.4935	062.4945
	3.5 mm	062.4936	062.4946
<b>Incl. Screw</b>		065.0036	

**Variobase® for Bridge/Bar Cylindrical, AH 3.5 mm**

		∅ 4.5 mm (RB/WB)
GH	1.5 mm	062.4961
<b>Incl. Screw</b>		065.0036

**Burn-out Copings for Variobase® for Crown**

		∅ 3.8 mm (RB/WB)	∅ 4.5 mm (RB/WB)
		065.0014	065.0015

**Burn-out Coping for Variobase® for Bridge/Bar Cylindrical**

		∅ 4.5 mm (RB/WB)
		065.0017

**Locator® Zest**

1 mmH	02326
2 mmH	02327
3 mmH	02328
4 mmH	02329
5 mmH	02330
6 mmH	02331

Up to 20° Divergence

Up to 40° Divergence


**Instruments and Accessories for Locator®**  
מפתחות ללוקי זסט

8519-10

8540-10

046.417/6

8393



**Plan Abutment for Screw-retained Abutments**

RB/WB Plan Abutment		POM	
	for Screw-retained Abutments, angulation 0°, gingiva height 1.5/2.5/3.5/4.5 mm		025.0073V4
	for Screw-retained Abutments, angulation 17°, gingiva height 3.5/4.5/5.5 mm		025.0074V4
	for Screw-retained Abutments, angulation 30°, gingiva height 3.5/4.5/5.5 mm		025.0075V4

#### Open-tray Closed-tray Analog Metal Scanbody Scanpost for Sirona® Cerec Digital Analog

16mm 065.0032

24mm 065.0034

13mm 065.4810

065.0022

065.0103

065.0038

065.0024

**Variobase® for Crown AS, AH 5.5 mm**

		∅ 4.5 mm (RB/WB)
GH	1.5 mm	062.4972
<b>Incl. Screw</b>		065.0037
	13 mm	046.786
	21 mm	046.787
	27 mm	046.788

**Burn-out Coping for Variobase® for Crown AS 25°**

		∅ 4.5 mm (RB/WB)
		065.0018

**Variobase® C for Sirona® Cerec**

		∅ 3.8 mm (RB/WB)	∅ 4.5 mm (RB/WB)
GH	1.5 mm	062.4981	062.4982
	<b>Incl. Screw</b>		065.0036

**Novaloc®**

		∅ 3.8 mm (RB/WB)
Angle	0°	—
	15°	—
GH	1.5 mm	062.4501
	2.5 mm	062.4502
	3.5 mm	062.4503
	4.5 mm	062.4504
	5.5 mm	062.4505
	6.5 mm	062.4506
7.5 mm	—	062.4512

**Processing Package titanium**

2010.601-NOV

**Novaloc®**

2010.731-NOV

2010.741-NOV

2010.751-NOV

**Processing Package PEEK**

2010.611-NOV

**SRA Screw-retained Abutments**

		∅ 4.6 mm (RB/WB)			
		Angle	0°	17°	30°
GH	1.5 mm		062.4722S		
	2.5 mm		062.4723S		
GH	3.5 mm		062.4724S	062.4733S	062.4743S
	4.5 mm		062.4725S	062.4734S	062.4744S
	5.5 mm		—	062.4735S	062.4745S

**SRA Burn-out Copings, Plastic**

		∅ 4.6 mm
Crown		023.4748
Bridge		023.4758
<b>Incl. Screw</b>		023.4763

**SRA Protective Caps**

		∅	4.6 mm
H	4.5 mm		024.0020-04
	5.1 mm		024.4323-04
H	6.6 mm		024.4324-04
	8.1 mm		024.4325-04

**Variobase® for Crown, AH 5.5 mm**

		∅ 5.5 mm (WB)
GH	0.75 mm	062.4953
	1.5 mm	062.4954
<b>Incl. Screw</b>		065.0036
<b>Burn-out Coping for Variobase® for Crown</b>		
		∅ 5.5 mm (WB)
		065.0016

**Variobase® for Bridge/Bar Cylindrical, AH 3.5 mm**

		∅ 4.5 mm (RB/WB)
GH	1.5 mm	062.4961
<b>Incl. Screw</b>		065.0036
<b>Burn-out Coping for Variobase® for Bridge/Bar Cylindrical</b>		
		∅ 4.5 mm (RB/WB)
		065.0017

**Anatomic Abutments**

		∅ 6.0 mm (RB/WB)		
		Angle	0°	17°
GH	2.5 mm		062.4103	062.4153
	3.5 mm		062.4104	062.4154
	<b>Incl. Screw</b>		065.0036	

**RB/WB Pre-milled Abutment Blank**

		∅
	11.5mm	Medentika Holder 062.4601
	15.8mm	Medentika Holder 062.4602
	12mm	Amann Girschbach 062.4603
<b>Incl. Screw</b>		065.0036

**Handpiece Screwdrivers**

	046.410
	046.411
	046.412

**Manual Screwdrivers**

046.400 13mm

046.401 21mm

046.402 27mm

**Ratchet and Torque Control Device**

046.119

046.049

**SRA Titanium Copings Regular**

		∅ 4.6 mm
Crown		023.4747
Bridge		023.4751
<b>Incl. Screw</b>		023.4763

**SRA Titanium Coping Slim**

		∅ 4.6 mm
Crown		024.0023
Bridge		024.0024
<b>Incl. Screw</b>		023.4763

**SRA Variobase® for Bridge/ Bar Cylindrical**

		∅ 4.6 mm
		023.0028
<b>incl. screw</b>		023.4763

**SRA Lab Processing Screw**

L 20mm	025.0006
L 10mm	025.0052

**SRA Occlusal Screw**

	023.4763
--	----------

**SRA Burn-out Coping for Variobase® for Bridge/Bar Cylindrical**

		∅ 4.6 mm
		023.0032

**Variobase® for Crown AS, AH 5.5 mm**

		∅ 5.5 mm (WB)
GH	1.5 mm	062.4971
<b>Incl. Screw</b>		065.0037
	13 mm	046.786
	21 mm	046.787
	27 mm	046.788

**Burn-out Coping for Variobase® for Crown AS 25°**

		∅ 5.5 mm (WB)
		065.0019

**WB Variobase® C for Sirona®**

		∅ 5.5 mm (WB)
GH	5.5 mm	062.4983
<b>Incl. Screw</b>		065.0036

**Abutment Basal + VB Basal Screw**

	065.0036
--	----------

**Basal Screw VB AS**

	065.0037
--	----------

#### Screw - Retained

**Abutment removal kit**

RB/WB Abutment-Removal Screw		for removal of damaged abutments, TAN	065.0007
Removal Tool		for BLX Basal Screw, left-hand, length 27 mm, stainless steel	065.0008
		for BLX Basal Screw, left-hand, length 21 mm, stainless steel	065.0009

#### Screw - Retained

**Abutment removal kit**

RB/WB Abutment-Removal Screw		for removal of damaged abutments, TAN	065.0007
Removal Tool		for BLX Basal Screw, left-hand, length 27 mm, stainless steel	065.0008
		for BLX Basal Screw, left-hand, length 21 mm, stainless steel	065.0009