

Straumann® Dental Implant System System Overview
BONE LEVEL PROSTHETICS
 MRK.032.Rev.04

Soft Tissue Management
 Temporary Restorations
 מבנים זמניים וכיפות ריפוי

Transfer Parts טרנספרים

Impression Posts
 Implant Analogs

טרנספרים ואנלוגים

Digital Impression Posts
 Scanbodies and Repositionable Analogs

פיני סריקה ואנלוגים דיגיטלים

Prosthetic Restoration חלקי שיקום

Titanium Abutments

מבנים מטילניום

Plastic Copings
 Titanium Copings

שרוליי יציקה
 בלנק טיטניום
 חלקי עק למולטי וינט (SRA)

Gold Abutments
 CoCr Abutments
 Locator
 לוקייטור, מבנה לייציקה
 קובלט כרום/זהב

Auxiliary Parts
 Screws

חלקים נלווים וברגים

SC Ø 2.9

Temporary Abutment for Crown

GH	Ø	Oval	
		1 mm	024.0011
2 mm	3 mm	024.0015	
		024.0016	

Implant-level Impression

Closed-tray Impression: 025.0020

Open-tray Impression: short 025.0021, long 025.0022

Analog: 025.0023

SC Metal Scanbody: 025.0077

SC Repositionable Implant Analog: 025.0024

NC

Temporary Abutments

Crown TAN	024.2371
Bridge TAN	024.2375
PMMA-Crown	024.2372

Implant-level Impression

Closed-tray Impression: 025.2201 L: 12 mm

Open-tray Impression: short 16 mm 025.2202, long 24 mm 025.2205

Analog: 025.2101

NC Metal Scanbody: 025.0078

Repositionable Implant Analog: 025.2102

Ø 3.3

Protective Caps

H	Ø	4.6 mm	
		5.1 mm	024.0020-04
5.1 mm	6.6 mm	024.4323-04	
		024.4324-04	
8.1 mm	024.4325-04		

Abutment-level Impression

Closed-tray Impression: Crown 025.2246, Bridge/Bar 025.0014

Open-tray Impression: Crown 025.2244, Bridge/Bar 025.0012

Analog: Ø 4.6 mm, < 0°, 023.4756

SRA Metal Scanbody: Ø 4.6 mm 025.0081

SRA Repositionable Implant Analog: 025.0008

RC

Temporary Abutments

Crown TAN	024.4371
Bridge TAN	024.4375
PMMA-Crown	024.4372

Implant-level Impression

Closed-tray Impression: 025.4201 L: 12 mm

Open-tray Impression: short 16 mm long 24 mm 025.4202, 025.4205

Analog: 025.4101

RC Metal Scanbody: 025.0079

Repositionable Implant Analog: 025.4102

RC

Protective Caps

H	Ø	4.6 mm	
		5.1 mm	024.0020-04
5.1 mm	6.6 mm	024.4323-04	
		024.4324-04	
8.1 mm	024.4325-04		

Abutment-level Impression

Closed-tray Impression: Crown 025.2246, Bridge/Bar 025.0014

Open-tray Impression: Crown 025.2244, Bridge/Bar 025.0012

Analog: Ø 4.6 mm, < 0°, 023.4756

SRA Metal Scanbody: Ø 4.6 mm 025.0081

SRA Repositionable Implant Analog: 025.0008

Kit BLT Modular קיט מודולרי

ST-BLT MODULAR

Screw-retained עבודה מוברגת

Variobase® for Crown

GH	Ø	Oval	
		1 mm	022.0038
2 mm	3 mm	022.0039	
		022.0040	

Cement-retained עבודה מודבקת

Anatomic Abutments

GH	Ø	0°		15°	
		2 mm	022.2102	022.2152	
3.5 mm	incl. screw	022.2104	022.2154	025.2900	

Standard Abutments

GH	Ø	0°		
		1.5	L 100	L 200-1
3.0	incl. screw	L 100-3	L 200-1-3	L 200-2-3

Screw-retained עבודה מוברגת

Variobase® for Crown

AH	3.5 mm	5.5 mm
GH 1 mm	025.2921	022.0027
GH 2 mm	022.0102	022.0106
GH 3 mm	022.0104	022.0108
incl. screw	025.2900	025.2900

Variobase® C Cerec Sirona (block S only)

GH 1 mm	022.0043
incl. screw	025.2900

Variobase® for Crown AS - Angulated Screw channel 25°

AH	3.5 mm	5.5 mm
GH	022.0084	022.0093
incl. screw	025.0053	025.0053

Variobase® for bridge/bar

AH 3.5 mm	022.0110
incl. screw	025.2926

Screw-retained עבודה מוברגת

Platform Ø 4.6

Screw Retained Abutment - (SRA) - עבודה מוברגת - מונים לעבודה מוברגת

GH	Ø	0°		17°		30°	
		1.5 mm	022.01285	-	-	-	-
2.5 mm	3.5 mm	Type A	022.01365	Type A	022.01445	Type A	022.01445
		Type B	022.01375	Type B	022.01465	Type B	022.01465
4.5 mm	5.5 mm	Type A	022.01405	Type A	022.01475	Type A	022.01485
		Type B	022.01415	Type B	022.01495	Type B	022.01495

Cement-retained עבודה מודבקת

Anatomic Abutments

GH	Ø	0°		15°	
		2 mm	022.4102	022.4152	022.4104
3.5 mm	incl. screw	022.4103	022.4107	025.4900	025.4900

Standard Abutments

GH	Ø	0°		
		1.5	L 110	L 210-1
3.0	incl. screw	L 110-3	L 210-1-3	L 210-2-3

Screw-retained עבודה מוברגת

Variobase® for Crown

AH	3.5 mm	5.5 mm
GH 1 mm	025.4921	022.0026
GH 2 mm	022.0103	022.0107
GH 3 mm	022.0105	022.0109
incl. screw	025.4900	025.4900

Variobase® C

GH 1 mm	022.0044
incl. screw	025.4900

Variobase® for Crown AS

AH	3.5 mm	5.5 mm
GH	022.0087	022.0096
incl. screw	025.0053	025.0053

Variobase® for bridge/bar

GH	022.0111
incl. screw	025.2926

Screw-retained עבודה מוברגת

Platform Ø 4.6

Screw Retained Abutment - (SRA) - עבודה מוברגת - מונים לעבודה מוברגת

GH	Ø	0°		17°		30°	
		1.5 mm	022.01325	-	-	-	-
2.5 mm	3.5 mm	Type A	022.01505	Type A	022.01585	Type A	022.01595
		Type B	022.01515	Type B	022.01615	Type B	022.01625
4.5 mm	5.5 mm	Type A	022.01545	Type A	022.01615	Type A	022.01625
		Type B	022.01555	Type B	022.01635	Type B	022.01635

NC Cementable Abutment CAD/CAM

Ø 3.5 mm	022.2311	NC Cementable Abut Ø 3.5/AH 4mm/GH 1mm
Ø 5.5 mm	022.2315	NC Cementable Abut Ø 3.5/AH 5.5mm/GH 1mm
Ø 4 mm	022.2321	NC Cementable Abut Ø 5/AH 4mm/GH 1mm
Ø 5 mm	022.2325	NC Cementable Abut Ø 5/AH 5.5mm/GH 1mm
Incl. Screw	025.2908	

Burn-out Copings for Variobase® for Crown

AH 3.5 mm	023.2756
AH 5.5 mm	023.0018

Burn-out Copings for Variobase® for Crown AS

Ø 4.6 mm	023.0025
----------	----------

Burn-out Copings for Variobase® for bridge/bar

Ø 4.6 mm	023.0029
----------	----------

SRA Titanium Coping

Ø 4.6 mm	023.4747
Crown	023.4751
Bridge	023.4758
incl. screw	024.0024

Variobase® Coping for Screw-retained Abutments

Ø 4.6 mm	023.0028
incl. screw	023.0032

RC Cementable Abutment CAD/CAM

Ø 4 mm	022.4321	RC Cementable Abut Ø 5/AH 4mm/GH 1mm
Ø 5 mm	022.4325	RC Cementable Abut Ø 5/AH 5.5mm/GH 1mm
Ø 6.5 mm	022.4331	RC Cementable Abut Ø 6.5/AH 4mm/GH 1mm
Ø 6.5 mm	022.4335	RC Cementable Abut Ø 6.5/AH 5.5mm/GH 1mm
Incl. Screw	025.4908	

Burn-out Copings for Variobase® for Crown

AH 3.5 mm	023.4759
AH 5.5 mm	023.0017

Burn-out Copings for Variobase® for Crown AS

Ø 4.6 mm	023.0026
----------	----------

Burn-out Copings for Variobase® for bridge/bar

Ø 4.6 mm	023.0030
----------	----------

SRA Titanium Coping

Ø 4.6 mm	023.4747
Crown	023.4751
Bridge	023.4758
incl. screw	024.0024

Variobase® Coping for Screw-retained Abutments

Ø 4.6 mm	023.0028
incl. screw	023.0032

RC Pre-milled Abutment Blank

Ø 11.5	Medentika Holder	022.0030
Ø 15.8	Medentika Holder	022.0031
Ø 12	Armann Girschbach Holder	022.0061

Gold Abutments

Crown	022.4410	incl. screw 025.4900
Bridge	022.4430	incl. screw 025.4902

CoCr

Crown	L 710
Bridge	L 710 R
incl. screw	L 60

SRA Titanium Coping

Ø 4.6 mm	023.4747
Crown	023.4748
Bridge	023.4758
incl. screw	024.0024

Variobase® Coping for Screw-retained Abutments

Ø 4.6 mm	023.0028
incl. screw	023.0032

Burn-out Copings for Variobase® for Screw-retained Abutments

Ø 4.6 mm	023.0032
----------	----------

Basal screw

Ø 2.9	025.0031
-------	----------

NC Pre-milled Abutment

Ø 11.5	Medentika Holder	022.0028
Ø 15.8	Medentika Holder	022.0029
Ø 12	Armann Girschbach Holder	022.0060

Gold Abutments

Crown	022.2410
incl. screw	025.2900

CoCr

Crown	L 700
Bridge	L 700 R
incl. screw	L 60

Gold Coping for SRA Ø 4.6

Ø 4.6 mm	023.4753
Crown	023.4754
Bridge	023.4754
incl. screw	023.4763

Occlusal Screw

Ø 3.3	023.4763
-------	----------

Lab Processing Screw

L 20mm	025.0006
L 10mm	025.0052

Locator® RC

1mmH	08170
2mmH	08171
3mmH	08172
4mmH	08173
5mmH	08174
6mmH	08175

Up to 20° Divergence: 8519-10

Up to 40° Divergence: 8540-10

Basal screw (VB+Gold)

VB Crown	025.4900
Bridge	025.4902

Basal screw (VB bridge)

Ø 2.9	025.2926
-------	----------

Basal screw AS

Ø 2.9	025.0053
-------	----------

Occlusal Screw for SRA Ø 4.6

Ø 3.3	023.4763
-------	----------

Lab Processing Screw for SRA Ø 4.6

L 20mm	025.0006
L 10mm	025.0052

Basal screw

Ø 2.9	025.0031
-------	----------

Locator® NC

1mmH	08422
2mmH	08423
3mmH	08424
4mmH	08425
5mmH	08425
6mmH	08426

Up to 20° Divergence: 8519-10

Up to 40° Divergence: 8540-10

Basal screw (VB+Gold)

Ø 2.9	025.2900
-------	----------

Basal screw (VB bridge)

Ø 2.9	025.2926
-------	----------

Basal screw AS (VB angulated)

Ø 2.9	025.0053
-------	----------

Occlusal Screw

Ø 3.3	023.4763
-------	----------

Lab Processing Screw

L 20mm	025.0006
L 10mm	025.0052

Locator® RC

1mmH	08170
2mmH	08171
3mmH	08172
4mmH	08173
5mmH	08174
6mmH	08175

Up to 20° Divergence: 8519-10

Up to 40° Divergence: 8540-10

Basal screw (VB+Gold)

VB Crown	025.4900
Bridge	025.4902

Basal screw (VB bridge)

Ø 2.9	025.2926
-------	----------

Basal screw AS

Ø 2.9	025.0053
-------	----------

Occlusal Screw for SRA Ø 4.6

Ø 3.3	023.4763
-------	----------

Lab Processing Screw for SRA Ø 4.6

L 20mm	025.0006
L 10mm	025.0052

Locator® לוקייטור

Instruments and Accessories for Locator®

נספחות לוקייטור

046.416

046.417

8393

Locator®

NC	RC	
	1mmH	08170
2mmH	08422	08171
3mmH	08423	08172
4mmH	08424	08173
5mmH	08425	08174
6mmH	08426	08175

Up to 20° Divergence: 8519-10

Up to 40° Divergence: 8540-10

15 lbs 3 lbs 5 lbs

20° 30°

0 lbs 1 lbs 2 lbs 4 lbs

Instruments and Accessories for Locator®

נספחות פרוטטיים

מכשיר עור לפרפת יוני

046.792

046.402 27 mm

046.401 21 mm

046.400 15 mm

046.119

046.049

נספחות פרוטטיים לזוויתן

046.412

046.411

046.410