

BONE LEVEL SURGERY

MRK.032.Rev.03

Ø 2.9 mm



**Bone Level Tapered (BLT)
SC**
Ø 2.9 mm



Titanium SLA®

6 mm	—
8 mm	—
10 mm	—
12 mm	—
14 mm	—
16 mm	—

Roxolid® SLActive®

6 mm	—
8 mm	—
10 mm	021.0010
12 mm	021.0012
14 mm	021.0014
16 mm	—

Round Burs and Needle Drills



Needle Drill
Ø 1.6 mm
026.0056

Pilot and Twist Drills



Pilot Drill 1
Ø 2.2 mm
066.1301 short
066.1701 long

**Alignment Pins
Depth Gauges**



Alignment Pin
Ø 2.2 mm
046.799

Profile Drills



BLT Profile Drill
Ø 2.9 mm
026.0061

Taps



BLT Tap for Adapter
Ø 2.9 mm
026.0063

**Closure Caps
Healing Abutments**

SC

Closure Caps



H 0.5 mm
024.00065

Healing Abutments



H 2 mm
024.00075



H 3.5 mm
024.00085



H 5 mm
024.00095

Ø 3.3 mm

Bone Level Tapered (BLT)

NC

Ø 3.3 mm



Ø 2.0 mm



6 mm

–

8 mm

021.3508

Roxolid®

10 mm

021.3510

12 mm

021.3512

SLA®

14 mm

021.3514

only

16 mm

021.3516

6 mm

–

8 mm

021.3308

10 mm

021.3310

12 mm

021.3312

14 mm

021.3314

16 mm

021.3316



Round Bur
Ø 2.3 mm
044.003



Round Bur
Ø 3.1 mm
044.004



Round Bur
Ø 2.3 mm
044.003



Pilot Drill
Ø 2.2 mm
066.1301 short
066.1701 long



X VeloDrill
Ø 2.8 mm
066.1302 short
066.1702 long



Pilot Drill
Ø 2.2 mm
066.1301 short
066.1701 long



Alignment Pin
Ø 2.2 mm
046.799



Depth Gauge
Ø 2.8 mm
046.800



Alignment Pin
Ø 2.2 mm
046.799



BLT Profile Drill
034.324



BLT Tap
034.345

NC

Closure Caps



H 0 mm
024.21005



Ø 3.6 mm
H 2 mm
024.22225



Ø 3.6 mm
H 3.5 mm
024.22245



Ø 3.6 mm
H 5 mm
024.22265



Ø 4.8 mm
H 2 mm
024.22425



Ø 4.8 mm
H 3.5 mm
024.22445

Healing Abutments

Ø 4.1 mm

Ø 4.8 mm

Bone Level Tapered (BLT)
RC
Ø 4.1 mm



—

021.5408

021.5410

021.5412

021.5414

021.5416

—

021.5308

021.5310

021.5312

021.5314

021.5316



Bone Level Tapered (BLT)
RC
Ø 4.8 mm



—

6 mm

8 mm

10 mm

12 mm

14 mm

16 mm

6 mm

8 mm

10 mm

12 mm

14 mm

16 mm

021.7408

021.7410

021.7412

021.7414

021.7416

—

021.7308

021.7310

021.7312

021.7314

021.7316



Round Bur
Ø 3.1 mm
044.004



Round Bur
Ø 2.3 mm
044.003



X VeloDrill
Ø 2.8 mm
066.1302 short
066.1702 long



X VeloDrill
Ø 3.5 mm
066.1304 short
066.1704 long



Pilot Drill
Ø 2.2 mm
066.1301 short
066.1701 long



X VeloDrill
Ø 2.8 mm
066.1302 short
066.1702 long



Depth Gauge
Ø 2.8 mm
046.800



Depth Gauge
Ø 3.5 mm
046.802



Alignment Pin
Ø 2.2 mm
046.799



Depth Gauge
Ø 2.8 mm
046.800



BLT Profile Drill
034.325



BLT Profile Drill
034.326



BLT Tap
034.346



BLT Tap
034.347

RC

RC

Closure Caps

Healing Abutments



H 0 mm
024.4100S



Ø 5.0 mm
H 2 mm
024.4222S






Ø 5.0 mm
H 4 mm
024.4224S



Ø 5.0 mm
H 6 mm
024.4226S



Ø 6 mm
H 6 mm
024.4228S



Basic Tray 041.778

Adapter for Ratchet 046.460 046.461

Adapter for Handpiece 046.470 046.471

Drill Extender 040.563

BLT Tap 026.0063

BLT Tap 034.345

BLT Tap 034.346

BLT Tap 034.347

BLT Profile Drill 026.0061

BLT Profile Drill 034.324

BLT Profile Drill 034.325

BLT Profile Drill 034.326

Round Bur 044.004 044.003

Needle Drill 026.0056

X VeloDrill™ 066.1301 066.1701

X VeloDrill™ 066.1302 066.1702

X VeloDrill™ 066.1304 066.1704

X VeloDrill™ 066.1306 066.1706

Alignment Pin 046.799

Depth Gauge 046.800

Depth Gauge 046.802

Depth Gauge 046.804

Ø 2.9, Ø 3.3, Ø 4.1, Ø 4.8

Ø 3.1, Ø 2.3, Ø 1.8

Ø 2.2, Ø 2.8, Ø 3.5, Ø 4.2

long, short

A Module 041.761

Ratchet Tray 041.766

Ratchet 046.119

Torque Control 046.049

Holding Key 046.064

Release Aid Tray 041.767

Release Aid N 026.2558

Release Aid R/W 026.4558

Release Aid S* 026.0073

SC Adapter* for Ratchet 025.0043

48 hour Tool SC* 026.0072

SCS Screwdriver for ratchet 046.400 046.401 046.402

SCS Screwdriver for Handpiece 046.410 046.411 046.412

BL Bone Profiler 026.0022 026.0023 026.0024

Bone Profiler Set Grommet Tray 6 small 041.762

48h tool Set Grommet Tray 6 large 041.763

48 hour Tool NC/RC 026.2048 026.4048

SC Adapter* for Handpiece 025.0042

Bone Profiler SC* 026.0070

SCS Screwdriver Set Grommet Tray 3 small + 3 large 041.764

BL Guiding Cylinder NC/RC 026.00255 026.00265

SC Auxiliary Set* Grommet Tray 3 small + 3 large 041.764