

Titanium SLA®

Roxolid® SLActive®

Round Burs and Needle Drills

Pilot and Twist Drills

**Alignment Pins
Depth Gauges**

Profile Drills

Taps

**Closure Caps
Healing Abutments**

Ø 2.9 mm

Bone Level Tapered (BLT)
SC
Ø 2.9 mm

Ø 1.8 mm

6 mm	-
8 mm	-
10 mm	-
12 mm	-
14 mm	-
16 mm	-
6 mm	-
8 mm	-
10 mm	021.0010
12 mm	021.0012
14 mm	021.0014
16 mm	-

Needle Drill
Ø 1.6 mm
026.0056

Pilot Drill 1
Ø 2.2 mm
066.1301 short
066.1701 long

Alignment Pin
Ø 2.2 mm
046.799

BLT Profile Drill
Ø 2.9 mm
026.0061

BLT Tap for Adapter
Ø 2.9 mm
026.0063

SC

Closure Caps

H 0.5 mm
024.00065

Healing Abutments

H 2 mm
024.00075

H 3.5 mm
024.00085

H 5 mm
024.00095

Ø 3.3 mm

Bone Level Tapered (BLT)
NC
Ø 3.3 mm

Ø 2.0 mm

6 mm	-
8 mm	021.3508
10 mm	021.3510
12 mm	021.3512
14 mm	021.3514
16 mm	021.3516
6 mm	-
8 mm	021.3308
10 mm	021.3310
12 mm	021.3312
14 mm	021.3314
16 mm	021.3316

Roxolid®
SLA®
only

Round Bur
Ø 2.3 mm
044.003

Round Bur
Ø 3.1 mm
044.004

Pilot Drill
Ø 2.2 mm
066.1301 short
066.1701 long

X VeloDrill
Ø 2.8 mm
066.1302 short
066.1702 long

Alignment Pin
Ø 2.2 mm
046.799

Depth Gauge
Ø 2.8 mm
046.800

BLT Profile Drill
034.324

BLT Tap
034.345

NC

Closure Caps

H 0 mm
024.21005

Healing Abutments

Ø 3.6 mm
H 2 mm
024.22225

Ø 3.6 mm
H 3.5 mm
024.22245

Ø 3.6 mm
H 5 mm
024.22265

Ø 4.8 mm
H 2 mm
024.22425

Ø 4.8 mm
H 3.5 mm
024.22445

Ø 4.1 mm

Bone Level Tapered (BLT)
RC
Ø 4.1 mm

Ø 2.7 mm

6 mm	-
8 mm	021.5408
10 mm	021.5410
12 mm	021.5412
14 mm	021.5414
16 mm	021.5416
6 mm	-
8 mm	021.5308
10 mm	021.5310
12 mm	021.5312
14 mm	021.5314
16 mm	021.5316

Round Bur
Ø 2.3 mm
044.003

Round Bur
Ø 3.1 mm
044.004

Pilot Drill
Ø 2.2 mm
066.1301 short
066.1701 long

X VeloDrill
Ø 2.8 mm
066.1302 short
066.1702 long

X VeloDrill
Ø 3.5 mm
066.1304 short
066.1704 long

Alignment Pin
Ø 2.2 mm
046.799

Depth Gauge
Ø 2.8 mm
046.800

Depth Gauge
Ø 3.5 mm
046.802

BLT Profile Drill
034.325

BLT Tap
034.346

RC

Closure Caps

H 0 mm
024.41005

Healing Abutments

Ø 5.0 mm
H 2 mm
024.42225

Ø 5.0 mm
H 4 mm
024.42245

Ø 5.0 mm
H 6 mm
024.42265

Ø 6.5 mm
H 2 mm
024.42425

Ø 6.5 mm
H 4 mm
024.42445

Ø 6.5 mm
H 6 mm
024.42465

Ø 4.8 mm

Bone Level Tapered (BLT)
RC
Ø 4.8 mm

Ø 3.2 mm

6 mm	-
8 mm	021.7408
10 mm	021.7410
12 mm	021.7412
14 mm	021.7414
16 mm	021.7416
6 mm	-
8 mm	021.7308
10 mm	021.7310
12 mm	021.7312
14 mm	021.7314
16 mm	021.7316

Round Bur
Ø 2.3 mm
044.003

Round Bur
Ø 3.1 mm
044.004

Pilot Drill
Ø 2.2 mm
066.1301 short
066.1701 long

X VeloDrill
Ø 2.8 mm
066.1302 short
066.1702 long

X VeloDrill
Ø 3.5 mm
066.1304 short
066.1704 long

X VeloDrill
Ø 4.2 mm
066.1306 short
066.1706 long

Alignment Pin
Ø 2.2 mm
046.799

Depth Gauge
Ø 2.8 mm
046.800

Depth Gauge
Ø 3.5 mm
046.802

Depth Gauge
Ø 4.2 mm
046.804

BLT Profile Drill
034.326

BLT Tap
034.347

RC

Closure Caps

Ø 5.0 mm
H 2 mm
024.42225

Ø 5.0 mm
H 4 mm
024.42245

Ø 5.0 mm
H 6 mm
024.42265

Ø 6.5 mm
H 2 mm
024.42425

Ø 6.5 mm
H 4 mm
024.42445

Ø 6.5 mm
H 6 mm
024.42465

Basic Tray 041.778

Adapter for Ratchet
046.460
046.461

Adapter for Handpiece
046.470
046.471

Drill Extender
040.563

BLT Profile Drill
026.0061

BLT Tap
034.345

BLT Tap
034.346

BLT Tap
034.347

BLT Profile Drill
034.324

BLT Profile Drill
034.325

BLT Profile Drill
034.326

Round Bur
044.004
044.003

Needle Drill
026.0056

X VeloDrill™
066.1301
066.1701

X VeloDrill™
066.1302
066.1702

X VeloDrill™
066.1304
066.1704

X VeloDrill™
066.1306
066.1706

Alignment Pin
046.799

Depth Gauge
046.800

Depth Gauge
046.802

Depth Gauge
046.804

A Module 041.761

Ratchet Tray
041.766

Release Aid Tray
041.767

Ratchet
046.719

Torque Control
046.049

Holding Key
046.064

Release Aid N
026.2558

Release Aid R/W
026.4558

Release Aid S*
026.0073

SCS Screwdriver for Ratchet
026.0022

BL Bone Profiler NC/RC
046.400
046.401
046.402

48 hour Tool Grommet Tray 6 small
046.410
046.411
046.412

BL Guiding Cylinder NC/RC
026.0025

48 hour Tool Grommet Tray 6 large
041.763

SC Adapter* for Handpiece Grommet Tray 3 small + 3 large
025.0043

48 hour Tool SC*
026.0072

SCS Screwdriver Set Grommet Tray 3 small + 3 large
041.764

BL Guiding Cylinder NC/RC
026.0025

SC Auxiliary Set* Grommet Tray 3 small + 3 large
041.764