

**Titanium SLA®**

**Roxolid® SLActive®**

**Round Burs and Needle Drills**

**Pilot and Twist Drills**

**Alignment Pins  
Depth Gauges**

**Profile Drills**

**Taps**

**Closure Caps  
Healing Abutments**

**Ø 2.9 mm**

Bone Level Tapered (BLT)  
SC  
Ø 2.9 mm  
Ø 1.8 mm

6 mm	-
8 mm	-
10 mm	-
12 mm	-
14 mm	-
16 mm	-
6 mm	-
8 mm	-
10 mm	021.0010
12 mm	021.0012
14 mm	021.0014
16 mm	-

Needle Drill  
Ø 1.6 mm  
026.0056

Pilot Drill 1  
Ø 2.2 mm  
066.1301 short  
066.1701 long

Alignment Pin  
Ø 2.2 mm  
046.799

BLT Profile Drill  
Ø 2.9 mm  
026.0061

BLT Tap for Adapter  
Ø 2.9 mm  
026.0063

**SC**

Closure Caps  
H 0.5 mm  
024.00065

Healing Abutments  
H 2 mm  
024.00075  
H 3.5 mm  
024.00085  
H 5 mm  
024.00095

**Ø 3.3 mm**

Bone Level Tapered (BLT)  
NC  
Ø 3.3 mm  
Ø 2.0 mm

6 mm	-
8 mm	021.3508
10 mm	021.3510
12 mm	021.3512
14 mm	021.3514
16 mm	021.3516
6 mm	-
8 mm	021.3308
10 mm	021.3310
12 mm	021.3312
14 mm	021.3314
16 mm	021.3316

Roxolid®  
SLA®  
only

Round Bur  
Ø 2.3 mm  
044.003

Round Bur  
Ø 3.1 mm  
044.004

Pilot Drill  
Ø 2.2 mm  
066.1301 short  
066.1701 long

X VeloDrill  
Ø 2.8 mm  
066.1302 short  
066.1702 long

Alignment Pin  
Ø 2.2 mm  
046.799

Depth Gauge  
Ø 2.8 mm  
046.800

BLT Profile Drill  
034.324

BLT Tap  
034.345

**NC**

Closure Caps  
H 0 mm  
024.21005

Healing Abutments  
Ø 3.6 mm  
H 2 mm  
024.22225  
Ø 3.6 mm  
H 3.5 mm  
024.22245  
Ø 3.6 mm  
H 5 mm  
024.22265  
Ø 4.8 mm  
H 2 mm  
024.22425  
Ø 4.8 mm  
H 3.5 mm  
024.22445

**Ø 4.1 mm**

Bone Level Tapered (BLT)  
RC  
Ø 4.1 mm  
Ø 2.7 mm

6 mm	-
8 mm	021.5408
10 mm	021.5410
12 mm	021.5412
14 mm	021.5414
16 mm	021.5416
6 mm	-
8 mm	021.5308
10 mm	021.5310
12 mm	021.5312
14 mm	021.5314
16 mm	021.5316

Round Bur  
Ø 2.3 mm  
044.003

Round Bur  
Ø 3.1 mm  
044.004

Pilot Drill  
Ø 2.2 mm  
066.1301 short  
066.1701 long

X VeloDrill  
Ø 2.8 mm  
066.1302 short  
066.1702 long

X VeloDrill  
Ø 3.5 mm  
066.1304 short  
066.1704 long

Alignment Pin  
Ø 2.2 mm  
046.799

Depth Gauge  
Ø 2.8 mm  
046.800

Depth Gauge  
Ø 3.5 mm  
046.802

BLT Profile Drill  
034.325

BLT Tap  
034.346

**RC**

Closure Caps  
H 0 mm  
024.41005

Healing Abutments  
Ø 5.0 mm  
H 2 mm  
024.42225  
Ø 5.0 mm  
H 4 mm  
024.42245  
Ø 5.0 mm  
H 6 mm  
024.42265  
Ø 6.5 mm  
H 2 mm  
024.42425  
Ø 6.5 mm  
H 4 mm  
024.42445  
Ø 6.5 mm  
H 6 mm  
024.42465

**Ø 4.8 mm**

Bone Level Tapered (BLT)  
RC  
Ø 4.8 mm  
Ø 3.2 mm

6 mm	-
8 mm	021.7408
10 mm	021.7410
12 mm	021.7412
14 mm	021.7414
16 mm	021.7416
6 mm	-
8 mm	021.7308
10 mm	021.7310
12 mm	021.7312
14 mm	021.7314
16 mm	021.7316

Round Bur  
Ø 2.3 mm  
044.003

Round Bur  
Ø 3.1 mm  
044.004

Pilot Drill  
Ø 2.2 mm  
066.1301 short  
066.1701 long

X VeloDrill  
Ø 2.8 mm  
066.1302 short  
066.1702 long

X VeloDrill  
Ø 3.5 mm  
066.1304 short  
066.1704 long

X VeloDrill  
Ø 4.2 mm  
066.1306 short  
066.1706 long

Alignment Pin  
Ø 2.2 mm  
046.799

Depth Gauge  
Ø 2.8 mm  
046.800

Depth Gauge  
Ø 3.5 mm  
046.802

Depth Gauge  
Ø 4.2 mm  
046.804

BLT Profile Drill  
034.326

BLT Tap  
034.347

**RC**

Closure Caps  
Ø 5.0 mm  
H 2 mm  
024.42225  
Ø 5.0 mm  
H 4 mm  
024.42245  
Ø 5.0 mm  
H 6 mm  
024.42265  
Ø 6.5 mm  
H 2 mm  
024.42425  
Ø 6.5 mm  
H 4 mm  
024.42445  
Ø 6.5 mm  
H 6 mm  
024.42465

**Basic Tray 041.778**

Adapter for Ratchet  
046.460  
046.461

Adapter for Handpiece  
046.470  
046.471

Drill Extender  
040.563

BLT Tap  
026.0063  
034.345  
034.346  
034.347

BLT Profile Drill  
034.324  
034.325  
034.326

Round Bur  
044.004  
044.003

Needle Drill  
026.0056

X VeloDrill™  
066.1301  
066.1701  
066.1302  
066.1702  
066.1304  
066.1704  
066.1306  
066.1706

**A Module 041.761**

Ratchet Tray  
041.766

Release Aid Tray  
041.767

Ratchet  
046.719

Torque Control  
046.049

Holding Key  
046.064

Release Aid N  
026.2558

Release Aid R/W  
026.4558

Release Aid S\*  
026.0073

SCS Screwdriver for Ratchet  
026.0022  
046.400  
046.401  
046.402

SCS Screwdriver for Handpiece  
046.410  
046.411  
046.412

Bone Profiler Set  
Grommet Tray 6 small  
041.762

BL Bone Profiler  
026.0023  
026.0024

48 hour Tool Set  
Grommet Tray 6 large  
041.763

48 hour Tool SC\*  
026.0072

SC Adapter\* for Handpiece  
026.0070

SC Adapter\* for Ratchet  
026.0043

BL Guiding Cylinder NC/RC  
026.0025  
041.764

SC Auxiliary Set\*  
Grommet Tray 3 small + 3 large  
041.764