




















































Straumann® Dental Implant System System Overview TISSUE LEVEL SURGERY MRK.034.Rev.03		Ø 3.3 mm				Ø 4.1 mm				Ø 4.8 mm																																																					
		Standard Plus (SP) NNC Ø 3.5 mm  Ø 3.3 mm		Standard Plus (SP) RN Ø 4.8 mm  Ø 3.3 mm		Standard Plus (SP) RN Ø 4.8 mm  Ø 4.1 mm		Standard Plus (SP) RN Ø 4.8 mm  Ø 4.8 mm		Standard Plus (SP) WN Ø 6.5 mm  Ø 4.8 mm																																																					
Titanium SLA®		4 mm	–	4 mm	–	4 mm	–	4 mm	–	4 mm	–	6 mm	–	6 mm	–	6 mm	–	6 mm	–	6 mm	–	8 mm	043.4165	Roxolid® SLA® only	8 mm	043.2605	Roxolid® SLA® only	8 mm	043.9265	8 mm	043.9405	8 mm	043.9495	10 mm	043.4175	10 mm	043.2615	10 mm	043.9275	10 mm	043.9415	10 mm	043.9505	12 mm	043.4185	12 mm	043.2625	12 mm	043.9285	12 mm	043.9425	12 mm	043.9515	14 mm	043.4195	14 mm	043.2635	14 mm	043.9295	14 mm	043.9435	14 mm	–
Roxolid® SLActive®		6 mm	–	6 mm	–	4 mm	033.0435	4 mm	033.0445	4 mm	033.0455	8 mm	033.4165	8 mm	033.5115	6 mm	033.5605	6 mm	033.5905	6 mm	033.6105	10 mm	033.4175		10 mm	033.5125		8 mm	033.5615	8 mm	033.5915	8 mm	033.6115	12 mm	033.4185	12 mm	033.5135	10 mm	033.5625	10 mm	033.5925	10 mm	033.6125	14 mm	033.4195	14 mm	033.5145	12 mm	033.5635	12 mm	033.5935	12 mm	033.6135	14 mm	–	14 mm	033.5645	14 mm	–	14 mm	–		
Round Burs		 Round Bur Ø 2.3 mm 044.003		→	 Round Bur Ø 3.1 mm 044.004		 Round Bur Ø 2.3 mm 044.003		→	 Round Bur Ø 3.1 mm 044.004		 Round Bur Ø 2.3 mm 044.003		→	 Round Bur Ø 3.1 mm 044.004																																																
Pilot and Twist Drills		 Pilot Drill 1 Ø 2.2 mm 044.210 short 044.211 long		→	 Pilot Drill 2 Ø 2.8 mm 044.214 short 044.215 long		 Pilot Drill 1 Ø 2.2 mm 044.210 short 044.211 long		→	 Pilot Drill 2 Ø 2.8 mm 044.214 short 044.215 long		→	 Twist Drill PRO Ø 3.5 mm 044.250 short 044.251 long		 Pilot Drill 1 Ø 2.2 mm 044.210 short 044.211 long		→	 Pilot Drill 2 Ø 2.8 mm 044.214 short 044.215 long		→	 Twist Drill PRO Ø 2.8 mm 044.250 short 044.251 long		→		 Twist Drill PRO Ø 4.2 mm 044.254 short 044.255 long																																						
Alignment Pins Depth Gauges		 Alignment Pin Ø 2.2 mm 046.799		→	 Depth Gauge Ø 22.8 mm 046.800		 Alignment Pin Ø 2.2 mm 046.799		→	 Depth Gauge Ø 2.8 mm 046.800		→	 Depth Gauge Ø 3.5 mm 046.802		 Alignment Pin Ø 2.2 mm 046.799		→	 Depth Gauge Ø 22.8 mm 046.800		→	 Depth Gauge Ø 3.5 mm 046.802		→		 Depth Gauge Ø 4.2 mm 046.804																																						
Profile Drills		 NNC Profile Drill 034.327		→	 SP Profile Drill 034.332		 SP Profile Drill 034.333		 SP Profile Drill 044.084																																																						
Taps		 NNC Tap Ø 3.3 mm 034.348		→	 SP Tap Ø 3.3 mm 034.351		 SP Tap Ø 4.1 mm 034.352		→	 BL Tap Ø 4.1 mm 034.349 רק לשתלי באורך 4 מ"מ		 SP Tap Ø 4.8 mm 034.353		 BL Tap Ø 4.1 mm 034.350 רק לשתלי באורך 4 מ"מ																																																	
Closure Caps Healing Caps		NNC		RN		RN		RN		WN																																																					
Closure Caps		Healing Abutments		Closure Cap - Flat		Healing Abutments		Closure Caps		Healing Abutments																																																					
 H 1.5 mm 048.3255		 H 3 mm 048.0715		 048.3715		 Ø 6 mm H 1.5 mm 048.3735		 048.3755		 H 2 mm 048.0385																																																					
 H 4.5 mm 048.0745				 Ø 5.5 mm H 2 mm 048.0335		 Ø 5.5 mm H 3 mm 048.0345		 H 3 mm 048.0395		 H 4.5 mm 048.0535																																																					
				 Ø 5.5 mm H 4.5 mm 048.0375																																																											